

Contamination Control Kit

ATP Test Kit



The Contamination Control Kit (CCK-4) is a high-precision test for measuring total and free ATP in liquid samples. It provides expanded flexibility over all-in-one test devices such as Aquasnap™ by allowing sample cleanup in the case of high quenching, particulates or fibers in aqueous samples. CCK-4 uses Hygiena's unique liquid-stable reagent which eliminates reconstitution of lyophilized pellets, multiple vials and extends the product shelf life. In addition CCK-4 comes with an ATP standard and 100 cuvettes making it a complete kit.

Key Benefits

- ✔ Detects microbial and nonmicrobial contamination
- ✔ Uses liquid-stable reagents versus lyophilized - 4 bottles versus 20 vials
- ✔ Extended shelf life - 12 months
- ✔ Offers flexibility when dealing with challenging samples
- ✔ Ideal solution for measuring contamination in process water, wastewater, water tanks and cooling towers
- ✔ Kit includes all components plus ATP standard
- ✔ Reduces cost by up to 30% and saves time

Additional Requirements

- ✔ Hygiena Luminometers
- ✔ Hold-Rite Applicator (HR002)
- ✔ Pipette and ATP-free tips for 100ul volume



Americas

941 Avenida Acaso
Camarillo, CA 93012
USA
Tel: 805.388.8007
Fax: 805.388.5531
Email: enquiries@hygiena.net
Web: www.hygiena.net

Europe

Unit 11 Wenta Business Centre
Colne Way Watford
Hertfordshire WD2 4ND UK
Tel: +44 (0)1923 818821
Fax: +44 (0)1923 818825
Email: enquiries@hygiena.net
Web: www.hygiena.net