

Snapshot™

Universal ATP Testing Device For All Charm Science Luminometers

Designed for companies conducting ATP monitoring with a Charm Science luminometer. Snapshot reduces costs and achieves more reproducible results with its unique patented design. Snapshot works the same as the test device it is replacing so there is no need to adjust Pass/Fail setting.



Key Features

- Simple, quick and user-friendly
- Unique liquid-stable reagent provides unprecedented accuracy and reproducibility
- All-in-one patented test design delivers more reproducible results by preserving entire sample
- 12-month shelf life at refrigerated temperatures (2-8°C)
- 4-week shelf life at room temperatures (21°C)
- Tolerant to temperature abuse
- No need to change procedures or Pass/Fail setting - Snapshot works the same as PocketSwab Plus ATP test
- Costs up to 40% less than Charm Pocketswab Plus ATP test

Snapshot works with the following luminometers:

- NovaLum™
- FireFly 1™
- FireFly 2™
- LUM-T™

**WORKS WITH
CHARM LUMINOMETERS AND
COSTS UP TO 40% LESS**



Europe

Unit 11 Wenta Business Centre
Colne Way Watford
Hertfordshire WD2 4ND UK
Tel: +44 (0)1923 818821
Fax: +44 (0)1923 818825
Email: enquiries@hygiene.net
Web: www.hygiene.net

Americas

941 Avenida Acaso
Camarillo, CA 93012
USA
Tel: 805.388.8007
Fax: 805.388.5531
Email: enquiries@hygiene.net
Web: www.hygiene.net

ORDER

NovaLum™	CH1616N
FireFly 2™	CH1616N
FireFly 1™	CH1616N
LUM-T™	CH1616N
Adapter Cap	CHCAP




Snapshot Instructions

Snapshot™


Part Number: CH1616N

For Use With: All Charm Luminometers


Cap sold separately (order# CHCAP)




Collect Sample




Re-Insert



Activate



Shake



Insert



Read

Looking For A Better System?



**TRADE-UP
PROGRAM**

Hygiene's SystemSURE Plus brings together new technologies with a durable design that gives users more sensitivity, reproducibility and robustness. In addition, it provides even more cost savings. To find out more about the SystemSURE Plus hygiene monitoring system contact Hygiene at enquiries@hygiene.net or call 888-HYGIENA.